


ALL200_200



Aluminum 200X200 mm Tk=15 mm 25X25 M6

Rating: Not Rated Yet

Price

Salesprice with discount

225,00 €

On Order

[Ask a question about this product](#)

Manufacturer [Linari Engineering srl](#)

Description

Our modular supports are necessary to prepare a safe, stable electrospinning workspace with precise and repeatable geometries over time.

These plates are made of aluminum and have precision perforations which permit the simple fastening of both the supports for the needles and the vertical supports for the flat collectors.

Thanks to the high electrical conductivity they can be used to create ground planes or high voltage distribution systems, for example for an array of coaxial needles.

All the components are interchangeable and can be connected together to create support structures for video cameras, pumps, lamps and customized collectors.

Dimensions

200X200 mm Tk=15 mm 25X25 M6

Units in box: 1