

COAX_L3M



Base for Modular coaxial needle - Layer 3: External

Rating: Not Rated Yet

Price

Salesprice with discount

860,00 €

7 Days

[Ask a question about this product](#)

Manufacturer [Linari Engineering srl](#)

Description

Base for modular coaxial needle layer 3 - External

Needles for electrospinning : COAX_L3M

The bases for the modular coaxial needles are made of AISI 316 stainless steel and form the support and the hydraulic connection for every specific layer of the system.

On this base can be installed every needles of family COAX-N3xx

To complete a coaxial system are need at least COAX_L1M and COAX_L2M with two needles as show below

BASE		NEEDLE			
code	Description	code	Ga	O.D	Pcs/ Conf
			uge	D	
COAX_L1M	Layer 1 - Internal	COA21	0.0	1.1	
		X_N	8	5	
		108	3	1	
		COA15	1.1	1.1	
		X_N	8	3	
		118	3	7	
COAX_L2M ^e	Layer 2 - Intermediate	COA21	0.0	1.1	
		X_N	8	5	
		208	3	1	
		COA15	1.1	1.1	
		X_N	8	3	
		218	3	7	
		COA11	3.2	1	
		X_N	0	3	
		230	5	9	
COAX_L3M	Layer 3 - external	COA15	1.1	1.1	
		X_N	8	3	
		318	3	7	
		COA11	3.2	1	
		X_N	0	3	
		330	5	9	

Needle not included

Units in box: 1