

COAX_N108



Tip for modular coaxial needle base: LAYER 1, GAUGE 21

Rating: Not Rated Yet

Price

Salesprice with discount

58,00 €

3-5 Days

[Ask a question about this product](#)

Manufacturer [Linari Engineering srl](#)

Description

The interchangeable tips make it easy to replace the tip of every layer of the coaxial needle in order to be able to adapt different diameters.

The assembly is done with a simple screwdriver, supplied as standard, thanks to a proprietary

connection of very high precision with the base, which makes it possible to achieve stable and repeatable geometries. All the tips are made of AISI 316 stainless steel and are laser welded for maximum cleanliness and precision.

Customized tips can be created upon request with diameters and lengths differing from those listed in this catalog.

Needle layer 1 O.D./I.D. 0.83/0.51mm GAUGE 21

To use with COAX_L1Mbase only.

Units in box: 1