

## COAX\_N308



Tip for modular coaxial needle base: Layer 3 Gauge 21

Rating: Not Rated Yet

**Price**

Salesprice with discount

58,00 €

3-5 Days

[Ask a question about this product](#)

Manufacturer [Linari Engineering srl](#)

### Description

The interchangeable tips make it easy to replace the tip of every layer of the coaxial needle in order to be able to adapt different diameters.

The assembly is done with a simple screwdriver, supplied as standard, thanks to a proprietary connection of very high precision with the base, which makes it possible to achieve stable and repeatable geometries. All the tips are made of AISI 316 stainless steel and are laser welded for maximum cleanliness and precision.

Customized tips can be created upon request with diameters and lengths differing from those listed in this catalog.

**Needles for electrospinning : COAX\_N308**

---

Needle layer 3 O.D./I.D.0.83/0.51 mm Gauge 21

To use with COAX\_L3M base only.

Units in box: 1