

## COAX\_SUP



Dielectric support for coaxial needles

Rating: Not Rated Yet

**Price**

Salesprice with discount

165,00 €

7 Days

[Ask a question about this product](#)

Manufacturer [Linari Engineering srl](#)

### Description

Supports in totally dielectric material allow the users to freely position their metal needles, connected to high voltage, inside the electrospinning chamber.

The cursor can scroll and rotate along the vertical bar. The needles, either two sealed coaxial needles (two or three layers) or a single modular coaxial needle, are fastened to the cursor.

There is also a version supporting up to two metal needles with Luer connector.

All the needles that are installed on the dielectric support must be connected to the respective syringes by flexible tubes. See the "tubes" section of this catalog.

Adjustable height from 50 to 200 mm from ground floor.

Materials: Black PVC and fiber glass.

Units in box: 1