



## Electrospinning Station



Electrospinning station hood

## Hood : Electrospinning Station

---

Rating: Not Rated Yet

### Price

Salesprice with discount

1600,00 €

*On Order*

[Ask a question about this product](#)

Manufacturer [Linari Engineering srl](#)

### Description

Esta is a compact electrospinning station, convenient and safe for operators and instrumentation, without risks of electric shock.

We improved the basic structure with an aluminum frame and a shelf to place generator board and control panel, to have all components needed by the electropinning station in a single place.

Our electrospinning hood is customized for chemistry use, with several accessories and special solutions design.

#### Safety devices

Door safety switch for HV	magnetic
Ground connection bar	inside hood
All metal part grounded	YES
Over pression relief trap	on top side
Safety window	5mm thk. plexiglass
Ventilation	40mm axial fan, until 10 m3/h

#### Dimensions & Features

Dimensions (LxWxH)	120 x 60 x 225 cm
Box material	Melamine wood
Overall height	120 cm
Internal dimension (LxWxH)	116 x 56x 60 cm

## Hood : Electrospinning Station

---

Flow rate control	YES
Window	Vertical sliding
Internal illumination	YES

### Available options:

- Polyethylene (PE) drip tray with 2(H) cm border to prevent liquid dispersion after spillage, white color
- ATEX fan 0,18 kW with variable speed controller (inverter) to discharge explosive vapours from hood

Units in box: 1