



Electrospinning Station



Electrospinning station hood

Hood : Electrospinning Station

Rating: Not Rated Yet

Price

Salesprice with discount

1600,00 €

On Order

[Ask a question about this product](#)

Manufacturer [Linari Engineering srl](#)

Description

Esta is a compact electrospinning station, convenient and safe for operators and instrumentation, without risks of electric shock.

We improved the basic structure with an aluminum frame and a shelf to place generator board and control panel, to have all components needed by the electrospinning station in a single place.

Our electrospinning hood is customized for chemistry use, with several accessories and special solutions design.

Safety devices

| | |
|---------------------------|-------------------------------|
| Door safety switch for HV | magnetic |
| Ground connection bar | inside hood |
| All metal part grounded | YES |
| Over pression relief trap | on top side |
| Safety window | 5mm thk. plexiglass |
| Ventilation | 40mm axial fan, until 10 m3/h |

Dimensions & Features

| | |
|----------------------------|-------------------|
| Dimensions (LxWxH) | 120 x 60 x 225 cm |
| Box material | Melamine wood |
| Overall height | 120 cm |
| Internal dimension (LxWxH) | 116 x 56x 60 cm |

Hood : Electrospinning Station

| | |
|-----------------------|------------------|
| Flow rate control | YES |
| Window | Vertical sliding |
| Internal illumination | YES |

Available options:

- Polyethylene (PE) drip tray with 2(H) cm border to prevent liquid dispersion after spillage, white color
- ATEX fan 0,18 kW with variable speed controller (inverter) to discharge explosive vapours from hood

Units in box: 1