

LLM-LLF-15



PVC tubing L=150 cm with Luer Lock M/F ends

Rating: Not Rated Yet

Price

Salesprice with discount

12,00 €

On Order

[Ask a question about this product](#)

Manufacturer [Linari Engineering srl](#)

Description

Either to use coaxial needles, both sealed and modular, or to create accurate needle/collector geometries, flexible tubes are always necessary to connect the needles to the relevant syringes.

We have developed some highly flexible, transparent PVC tubes with male and female Luer Lock connectors on the ends to ease this connection.

The various lengths available are ideal to create solution feeding extensions whether the

needles are fastened or moving.

Description

Pipe Material	PVC
Luer Lock connector s- type	Male/Female
OD/ID	4,1/3,0 mm
Length	1500 mm

Units in box: 1