

## NEL\_1020\_12PC



Luer needles - OD=1.0 mm L=20 mm flat - 12 pcs

Rating: Not Rated Yet

**Price**

Salesprice with discount

25,00 €

3-5 Days

[Ask a question about this product](#)

Manufacturer [Socorex Isba S.A.](#)

Description

The AISI 316 stainless steel needles with Luer connection are an essential accessory for any

electrospinning activity.

The flat tip needles, compared to those with bevel tip, permit a more accurate control of the electrical field guaranteeing its absolute symmetry and reducing the local ionization (Corona effect) of the tip which can alter the solution when using high voltage values.

**Description**

Connector material	Metal
Connector type	Luer lock - Female
Needle material	stainless steel
Needle OD	1.0 mm
Needle length	20 mm
Tip type	Flat

Units in box: 1