


PVC300_600



PVC 300X600 mm Tk=15 mm 25X25 M6

Rating: Not Rated Yet

Price

Salesprice with discount

784,00 €

On Order

[Ask a question about this product](#)

Manufacturer [Linari Engineering srl](#)

Description

Our modular supports are necessary to prepare a safe, stable workspace with precise and repeatable geometries over time.

These plates are made of PVC and have precision perforations which permit the simple fastening of both the supports for the needles and the vertical supports for the flat collectors.

They are the most widely used supports due to their good quality/price ratio and are an ideal solution in most applications.

All the components are interchangeable and can be connected together to create support structures for video cameras, pumps, lamps and customized collectors.

Dimensions

300X600 mm Tk=15 mm 25X25 M6

Units in box: 1