

SALL100_150



Aluminum 90° SUPPORT 100X150 mm Tk=3 mm

Rating: Not Rated Yet

Price

Salesprice with discount

35,00 €

On Order

[Ask a question about this product](#)

Manufacturer [Linari Engineering srl](#)

Description

The vertical supports enable the rigid or flexible flat collectors to be fastened into position in front of the needle for every static electrospinning application.

Made of aluminum, they feature a perforation allowing to install them on the modular supports with special screws for a secure and stable hold that guarantees electrospinning with accurate and repeatable geometries over time.

The aluminum supports make the electric field uniform and permit using also dielectric collectors such as films and polymeric networks.

Dimensions (vertical face and thickness)

100X150 mm Tk=3 mm

Units in box: 1