

## SLEX100\_150



Plexiglas 90° SUPPORT 100X150 mm Tk=5 mm

Rating: Not Rated Yet

**Price**

Salesprice with discount

30,00 €

*On Order*

[Ask a question about this product](#)

Manufacturer [Linari Engineering srl](#)

### Description

The vertical supports enable the rigid or flexible flat collectors to be fastened into position in front of the needle for every static electrospinning application.

Made of Plexiglas, they feature a perforation allowing to install them on the modular supports with special screws for a secure and stable hold that guarantees electrospinning with accurate

and repeatable geometries over time.

The use of dielectric supports in Plexiglas allows to correctly manage electricity-conducting collectors -with even small shapes and geometries- where to concentrate the electric field, and so the nanofibers, avoiding their being attracted by unwanted metal surfaces.

**Dimensions (vertical face and thickness)**

100X150 mm Tk=5 mm

Units in box: 1