

## SSCY\_80X120



Aluminum, thin wall drum D=80, Length=120 mm

Rating: Not Rated Yet

**Price**

Salesprice with discount

286,00 €

3-5 Days

[Ask a question about this product](#)

Manufacturer [Linari Engineering srl](#)

### Description

The drum collectors are made of high resistance hollow aluminum to reduce to a minimum both the weight and the moment of inertia so that they can be used at high speeds of rotation without triggering dangerous vibrations.

This collector and all the ones with length up to 150 mm can be used on all the rotation systems of Linari Engineering.

Every drum collector comes with bare, non-anodized and untreated aluminum outer surfaces to reduce any possible electrical resistance which would alter the electrical field making the electrospinning process non-reproducible.

**Dimensions**

Outer diameter	80 mm
Grip diameter	3.0 mm
Lenght	120 mm

**Packaging**

1 pcs/conf

Units in box: 1