

## SSWCY\_80X120



Stainless Steel wire drum D=80, Length=120 mm, N.12 wires with diameter=3 mm

Rating: Not Rated Yet

**Price**

Salesprice with discount

550,00 €

1-4 Weeks

[Ask a question about this product](#)

Manufacturer [Linari Engineering srl](#)

### Description

The drum collectors is made of high resistance stainless steel to increase mechanical resistance to be used at high speeds of rotation without triggering dangerous vibrations.

This collector can be used only with EasyDrum rotation systems by Linari Engineering.

This kind of collector is suitable to produce single aligned fibers and to create tissues with few contact point with collector for an easy handling during detachment.

**Dimensions**

Outer wire diameter	80 mm
Grip diameter	10.0 mm
Wire length	120 mm
Wire diameter	3.0 mm
Wire number	12

Units in box: 1